

[Home Page](#)[Cyber Specials](#)[What's New](#)[2001 catalog](#)[Instrumentation](#)[Suppliers](#)[How to Order](#)[E-mail Us](#)

Carbowax 20M Non-Bonded Phases

- Polyethylene glycol, MW 20,000
- Equivalent to USP Phase G16

Unlike DB-WAX, CAM, Carbowax 20M and OV-351 GC columns are not bonded and cross-linked, therefore solvent rinsing is not recommended .

Test Mix:

Part Number: 200-0020

Similar Phases:

DB-WAX, HP20M

CARBOWAX 20M ORDER GUIDE:

Inner Diameter (mm I.D.)	Length (meter)	Film Thickness (µm)	Temperature Limits (°C)	Part Number
0.20	25	0.10	60 to 220/240	118-2021
0.20	50	0.10	60 to 220/240	118-2051
0.25	15	0.25	60 to 220/240	112-2012
0.25	30	0.25	60 to 220/240	112-2032
0.25	60	0.25	60 to 220/240	112-2062
0.32	15	0.25	60 to 220/240	113-2012
0.32	25	0.30	60 to 220/240	113-2023
0.32	30	0.25	60 to 220/240	113-2032
0.32	50	0.30	60 to 220/240	113-2053
0.32	60	0.25	60 to 220/240	113-2062
0.53	30	0.50	60 to 200/240	115-2037

[Click here to return to the Fused Silica Product Page](#)

your chromatography specialist